

2 in 1 BODY THERAPY

User Manual



NEW CONCEPT OF KEEP FIT & SCULPTURE BODY

- Relieve shoulder stress
- Beautify back
- Sculpture waist line
- Raise the buttocks
- Beautify leg

This product won patents in many countries, and cooperated with the Ministry of Economic Affairs.

CONTENTS




NAME	2 in 1 BODY THERAPY
ITEM	CH-022A
POWER	AC100V-240V DC12V/1.5A
WATT	9.6W
CURRENT	0.8A
WEIGHT	280g
MATERIAL	ABS+Silica
ORIGIN	MADE IN TAIWAN

- ① Host Body + 38 CM Belt x 1 SET
- ② Host Body + 48 CM Belt x 1 SET
- ③ Adapter x 1 PC ④ Cable x 1 PC
- ⑤ Bottle X 1 PC (optional)

EMS FUNCTION


Function	Frequency	Instruction	EMS Intensity Instruction
P1	HI	5 Plus/Cycle	
P2	Low	1 Plus/Cycle	
P3	Middle	3 Plus/Cycle	
P4	Combination		
P5	Combination		
P6	Continuous		

Single Vibration Belt Connected




After connected belt, you can adjust the belt to fit your size.

Handheld-Vibration Belt Connected



After follow the step of left picture, turn the belt to the other side to be handheld vibration.

Double Vibration Connected




Operation Method


1. ON/OFF : 

Push the ON/OFF button with light in the bottom. The product will automatically perform computer setup speed function 1~5. Speed functions change in one minute and automatically shut down after 15 minutes.





2. Hi/Low :  

Push the Hi/Low button to adjust the speed of massage, you can choose the speed you want. The speed is separated to five grades.


3. EMS Function : 

User could start massage function at the same time demand for EMS activated. User push EMS-Hi button  to start EMS function with light in the bottom. The product will automatically perform computer setup EMS function P1. User could push EMS Function to change the functions. The EMS function is separated to six functions.

4. EMS Hi/Low :  

EMS speed is separated to S1 to S10. The first push  is to start EMS functions. Then, user could change the speed to powerful speed by push EMS Hi . On the contrary, user could push EMS Low  to change speed to low. And the last push  is to close the function of EMS. Every EMS functions you change, computer will automatically change speed into S1.

5. The product will shut down automatically if the use of time up to 15 minutes or the higher of the temperature. When starting EMS function, it will detect whether the product contact with the skin within 20 seconds. If there is no connection or bed connection, it will shut down automatically.



AKEE 【2 In 1 BODY THERAPY】 According to the needs of various parts of your body sculpting, It use the power of Hula hoop as the centrifugal powerful rotary massage which could synchronized with the digestive system. It help body gain natural detoxification and purification. And join the EMS (electrical muscle stimulation) function which allowing the body to gain the best motion effects.

AKEE 【2 In 1 BODY THERAPY】 It provides professional massage function as spa to break up excess fat.

AKEE 【2 In 1 BODY THERAPY】

- 1.Exclusive use of single motor patent structure roundabout design, so that the volume of product is smaller which help users place in any body parts and reduce the large burden of users.
- 2.The powerful rotary vibration massage is the best to keep fit body.
- 3.Computer multi-function control which makes products more easy to use and ergonomic.
- 4.It could be use in different body parts, also could be use as handheld massager
- 5.For the use of a large area, just connect two groups of massage belt to increase the scope of body parts.
- 6.The electronic transformer is portable to all overthe world.
- 7.This product is produced by industry, government and academia collaboration to establish the best commercialization.
- 8.Products have lightweight and significant effect.
- 9.Combination of EMS and massage functions help users achieve the movement and body sculpture effectively.

Safety and Precautions

- 1.Do not use within one hour before and after the meal, after use please at least drink 350cc warm water. It promote metabolism.
- 2.If feel discomfort, please stop using.
- 3.After using,please wipe clean EMS surface with cotton cloth. To avoid damaging, please avoid direct water rinse or wipe with a poor quality cloth.
- 4.EMS surfaces must be well contact with the skin. Because, poor contact will cause too much impedance differences. Before using, you can use the container with warter to spray on massage part.
- 5.If health concerns arise while using this product, please stop use and contact dealer.

AKEE [2 in 1 BODY THERAPY] Recommend the use of parts

Recommend the use parts of a Single Massage Belt



Arm

- Sculpture arms
Slim flabby



Thigh

- Eliminate fatigue
Sculpture thigh
Slender lines of thigh



Leg

- Promote leg circulation &
Eliminate fatigue



Body

- As hand-massager use in every part of body

Recommend the use parts of Double Massage Belt with connecting cable



Waist

- Pat and rotary massage &
Slender perfect curve



Abdomen

- Tight muscles
Break up excess fat



Thigh

- Eliminate fatigue
Sculpture thigh
Slender lines of thigh



Hip

- To improve
Qiaotun curve &
Tight muscles round

Handheld Massager



It's used for body massage to eliminate fatigue.



Thank you for purchasing AKEE 【2 in 1 BODY HERAPY】
Before you use 【2 in 1 BODY THERAPY】 , please read manual clearly. You need to fully understand the correct way to use, especially pay attention to the part of the "Safety and Precautions" Please keep user manual for future reference use.

ATTENTION

- ✘ The product is recommended to start from low speed, and then gradually increase the speed.
- ✘ The product is recommended to the use of time up to 15 minutes in any body parts. It avoid excessive use which can cause muscle soreness.
- ✘ If you feel itchy after use, it's not side effect. Because the body of the vibration generate physiological compensation effect, please ease to use.
- ✘ Do not use the product more than 30 minutes to maintain the life of product each time.
- ✘ Before and after use, to promote metabolism, please drink a large glass of water (300cc-500cc).
- ✘ If user has specific diseases, please consult a professional physician prior to use. Users of pregnancy, children and pace maker do not use the product.
- ✘ If health concerns arise while using this product, please stop use and contact seller.

Warranty Card

Sincerely thank you for your patronizing. This product is subject to stringent quality control. You can safely use ◦ In order to ensure product quality and protect your consumer rights, we will provide one year of free service guarantee for you.

Warranty shall be the host body to enjoy one year of free ✘ service from date of purchase.

The card if no detailed dealer stamp of the date of ✘ purchase and the Company stamp, as the invalid.

In the following circumstances, the Company will charge ✘ the materials fee.

1. Over the guarantee period.
2. Did not show the Warranty Card.
3. The contents of Warranty Card with the present matter is different.

It will not be replaced the warranty card, if you lost, please ✘ save properly.

Product Name AKEE 2 in 1 BODY THERAPY

Date _____

Dealer _____

Dealer's Signature

